

## **Freeze Dryer**



Advanced Materials Research Centre (AMRC), Kamand Campus Indian Institute of Technology Mandi, Himachal Pradesh-175075 Email: amrcoffice@iitmandi.ac.in

## **Only For Internal User**

□ User:

Date:

- □ Supervisor:
- □ Type of vials required: 1.) 40ml ( ) 2.) 750ml ( )
- □ No. of vails required:
- □ Name of the solvent in your samples.:

□ How many hours have your samples been frozen at -80 degrees.:

- □ For how may hours do you need slot: (12hr/24hr):
- □ Name of the lyophilizer want to use:
  1) 4.5L (-105 C)
  2.) 2.5L (-84 C)

## **Note:**

- Please make sure your samples are properly frozen at -80 degrees for at least 12 hours before using the instrument,
- After the slot is complete, turn off the instrument & clean it properly, make it dry and return the vials to the AMRC office as soon as you finish your slot.
- Don't use the DMSO in any of the lyophilizer.
- $\hfill\square$  In the 2.5L lyophilizer , only water is used as solvent.